

VISION Pro Automatic ESR Analyzer



Innovative Technology

- Direct testing with EDTA tubes
- Quick infrared scanning
- Good correlation with Westergren method
- Results in 20 minutes

Advanced Features

- LCD touch screen
- Dynamic sedimentation curve available
- Internal barcode reader and Thermal printer
- Flexible test times allocation via master and secondary dongle

Patented in U.S.



Specifications

Test Speed	20 min/cycle
Model Types	8 /16/32 sample positions
Dongle Types	Master dongle and secondary dongle
Test Tube Requirement	EDTA tube with a diameter of 12/13mm
Sample Volume	1.5ml-3.5ml
Display	7 inches LCD touch screen
Barcode Reader	Internal barcode reader
Printer	Built-in Thermal printer
LIS	Bidirectional connect
ESR Dynamic Curve	Yes
Internal QC Management	Yes
Temperature Correction	Yes
Internal Mixing	Yes
Power Supply	100-240V, 50/60Hz
Input Power Supply	220VA
Dimension	352mm×365mm×325mm
Weight	10.5kg

DISTRIBUTOR:

P/N: ENG-ENG-C009-1801



YHLO

SHENZHEN YHLO BIOTECH CO., LTD.

Headquarter: 13th Floor, Building A1, Kexing Science Park, No.15 Keyuan Rd., Nanshan District, Shenzhen 518057, P.R.China
Factory: No.5 Building, Lishan Industrial Area, Xinghai Road, Nanshan District, Shenzhen 518054, P.R.China

Tel: +86 755 2647 3359 - 8151
Fax: +86 755 2647 3319
Email: marketing@szyhlo.com
Website: www.szyhlo.com